

DiaLog

A Diabetic Tracking Solution

Team Members

-Eduardo Reyes | Project Manager, Front End

Braden Post | Back End Engineer

Aaron Massey | Front End Engineer

Andres Ventura | Back End Engineer

Saul Cubillo | Back End Engineer | Quality Assurance

Purpose

The ideal is to create an easy to use tool, one that those with diabetes can easily input their data, and get an accurate reading to their own body.

Since blood sugar varies depending on the users body, if they've eaten recently, how they've done along the week. We wanted to make a dynamic tool that can give results depending on each individual user.

Project Timeline

- Week 1: Select groups and identify project goals
- Week 2: Identify project scope, including target audience, format, and platform to be used
- Week 3: Finalize on roles and languages to be used; begin researching software requirements for completing the project
- Week 4: Develop use case diagrams and finalize on class models
- Week 5: Create sequence and data flow diagrams
- Week 7: Finalize on the early iterations of documentation
- Week 8: Program login logic and its interaction with user database
- Week 9: Create user interface and optimize user experience
- Week 10: Update documentation to reflect changes made to software
- Week 11: Debug application and present demonstration

Ease of Use



We designed our app to incorporate a one handed use case. Hence it was made with easy to reach buttons in the interface, so the user can quickly grab their phone after taking a blood reading, and immediately get their calculations and report on how their week has been, and if any abnormalities have been recorded.

Language I Swift


The language we chose for this project is swift for the following reasons:

- The simplicity of integrating with apple devices for ease of use
- Swifts many in-built functions and extensions to enhance development
- The ease of testing through X-code
- It compliments our project with a clean UI and navigation experience


Diagrams




User Interface | Login



The image shows a mobile app interface for a login screen. At the top, the word "DiaLog" is displayed in a red, sans-serif font. Below the title, there are two input fields: "Email" and "Password", both with light gray borders. Underneath the input fields, there are two buttons: "Sign Up" and "Log In", both with a light gray background and rounded corners. The entire interface is set against a white background within a black smartphone frame.



The image shows the same mobile app interface as the first one, but with a modal dialog box displayed in the center. The dialog box has a white background and a thin gray border. It contains the text "LOGIN ERROR: The email address is already in use by another account." in a black, sans-serif font. Below the text, there is a single button labeled "OK" in blue. The background of the app is dimmed to a light gray color. The "Sign Up" and "Log In" buttons are now red with white text. The "Email" input field contains the text "test1".




The image shows the same mobile app interface as the first one, but with a modal dialog box displayed in the center. The dialog box has a white background and a thin gray border. It contains the text "LOGIN ERROR: The password is invalid or the user does not have a password." in a black, sans-serif font. Below the text, there is a single button labeled "OK" in blue. The background of the app is dimmed to a light gray color. The "Sign Up" and "Log In" buttons are now red with white text. The "Email" input field contains the text "test1".


Login Functionality | Firebase

Firebase is a backend database that is specialized for mobile app development. It allows for easy integration between Unity, Web, IOS, and Android applications.

Allows for user authentication using email, phone number, Facebook and Google accounts.



Home Screen



Sign Out

Calculator

Have you eaten in the last 8 hours?

☐ Yes ☐ No


What is your current mmol/L reading? :

mmol/L

Calculate

Person icon

Home indicator



Sign Out

Calculator

Have you eaten in the last 8 hours?

☒ Yes ☐ No

What is your current mmol/L reading? :

4.3

Calculate

Your current Blood Glucose is (mg/dL):
77.39... mg/dL

Your blood sugar is normal

Person icon

Home indicator

The main screen begins at the calculator for ease of use, incase the user only needs the reading.

Then the tabs along the bottom of the screen can navigate to the user data, and the weekly report and graph.

Blood Glucose Calculations

The math mostly revolves around the mmol/L to mg/dL

```
print("Please enter mmol/L reading")
var num1 = Float(readLine())!!

var bloodglucose = mmol * 18
print("Your current blood glucose is \$(sum) mg/dL")
```

The calculation is also prepared to ask if the user has eaten recently, as that varies the output greatly.

```
if let inputData = readLine() {
    if inputData.lowercased() == "yes"{
        💡 print("subject is fasting")
    }
}
```

Benefits of mg/dL over mmol

There is great benefit to doing these calculations instead of taking the Blood glucose monitors mmol reading at face value.


Mmol/L provides the molarity (the number of molecules of any substance within a defined volume mostly used for fluids), while, mg/dL provides the concentration by the ratio of weight to volume, for example, milligrams per deciliter.

Mg/dL is the actual amount of blood glucose in the blood per liter, which is a more accurate calculation when factoring in the users health, biometrics, and diet.

Graphing Capability

The application will have the capability to display a user's blood sugar levels over the course of the past week in order to observe trends


More functionality is planned in order to adjust visibility of time periods and provide comments and insights on current trends



Biometric Updating (UX)

Upon creating an account, users will be prompted to provide their relevant biometric information

This information can be updated when necessary to maintain accurate calculation values



The image shows a mobile application screen with a white background and rounded corners. At the top left, there is a blue back arrow and the text '< Back'. In the center, the text 'Your Information:' is displayed. Below this, there is a vertical stack of six text input fields with light gray borders. The labels for these fields are 'First Name', 'Last Name', 'Weight', 'Height', 'Diabetes Type', and 'BMI'. At the bottom center of the form, there is a red rectangular button with the white text 'Save'.

